





Backing:1/8" white, crosslinked polyethylene
foamAdhesive:Non-irritating acrylic (M111)Liner:2 mil clear BOPP film, silicone treated on

Applications

- Electrodes: ECG, EKG, TENS
- AED defibrillator pads

Medical device mounting

one side

• Transdermal delivery system

Features	Benefits
Polyethylene Foam Backing	 Soft and flexible for maximum patient comfort Conformable and easily adapts to medical devices and body contours
Acrylic Adhesive	 Strong adhesion to surfaces
Polypropylene Liner	Excellent lay flat and dimensional stabilityEasily die-cut

Typical Test Values

	Standard UOM	Metric UOM	Test Method
Total Thickness w/liner	129 mils	3.277 mm	PSTC 133
Adhesive Thickness	2.0 mils	0.051 mm	PSTC 133
Carrier Thickness	125.0 mils	3.175 mm	PSTC 133
Liner Thickness	2.0 mils	0.051 mm	PSTC 133
Adhesion to Steel	110 oz/in	12.03 N/cm	PSTC 101 Mod
	Produc	t Safety	

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be noncytotoxic, and non-irritating.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60 °F - 80 °F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY assumes any legal responsibility.



25 Forge Pkwy Franklin, MA 02038 1-800-343-7875 Fax 1-800-328-4822 www.berryglobal.com