





Backing:1/32" white, crosslinked polyethylene foamAdhesive:Non-Sensitizing Acrylic (M108)Liner:74lb bleached polycoated kraft, silicone
treated on one side

- Applications
- Electrodes: ECG, EKG, TENS
- AED defibrillator pads

- Medical device mounting
- Transdermal delivery system

Features	Benefits
Polyethylene Foam Backing	 Soft and flexible for maximum patient comfort Conformable and easily adapts to medical devices and body contours
Acrylic Adhesive	High tack bonds well to various materialsStable to EtO and gamma sterilization
• 74lb Kraft Liner	Easily die-cut into various shapes and sizesGood dimensional stability and lay flat

Typical Test Values

		•					
	Standard UOM		Metric UOM		Test Method		
Total Thickness w/liner	38.8	mils	0.986	mm	PSTC 133		
Adhesive Thickness	2.0	mils	0.051	mm	PSTC 133		
Carrier Thickness	31.2	mils	0.792	mm	PSTC 133		
Liner Thickness	5.6	mils	0.142	mm	PSTC 133		
Adhesion to Steel	75	oz/in	8.21	N/cm	PSTC 101 Mod		
Product Safety							

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be noncytotoxic, non-irritating, and non-sensitizing.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60°F-80°F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY assumes any legal responsibility.



25 Forge Pkwy Franklin, MA 02038 1-800-343-7875 Fax 1-800-328-4822 www.berryglobal.com