



	Release Liner					
	, Adhesive	Adhesive:	Medical-grade acrylic (M118)			
		Liner:	80lb white kraft liner, silicone release treated on both sides with differential release			
Applications						
Medical device attachmeMedical device assembly		 Polyk 	bag closure			

Features	Benefits
Medical Grade Acrylic	 Good humidity and stress resistance Aggressive bond
	 Adheres well to foams and plastics including PE, TPO, ABS, and PP
80lb Release Liner	Easy liner removal

Typical Test Values

	Standard UOM	Metric UOM	Test Method				
Total Thickness w/liner	6.5 mils	0.165 mm	PSTC 133				
Adhesive Thickness	2.0 mils	0.051 mm	PSTC 133				
Liner Thickness	4.5 mils	0.114 mm	PSTC 133				
Adhesion to Steel	100 oz/in	10.94 N/cm	PSTC 101 Mod				
Product Safety							

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be noncytotoxic and non-irritating.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60°F-80°F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY assumes any legal responsibility.



25 Forge Pkwy Franklin, MA 02038 1-800-343-7875 Fax 1-800-328-4822 www.berryglobal.com