



Release Liner			
, Adhesive	Adhesive:	Medical-grade acrylic (M116)	
	Liner:	74lb bleached polycoated kraft, silicone trated on both sides	
	Applications		

Medical device attachment

• Medical device assembly

Features	Benefits		
Acrylic Adhesive	 Excellent adhesion to skin 		
	 Good balance of shear, tack, and adhesion 		
74lb Release Liner	 Good dimensional stability and lay-flat 		
	Easy to die-cut		

Typical Test Values

	Standard UOM	Metric UOM	Test Method
Total Thickness w/liner	7.0 mils	0.178 mm	PSTC 133
Adhesive Thickness	1.1 mils	0.028 mm	PSTC 133
Liner Thickness	5.9 mils	0.150 mm	PSTC 133
Adhesion to Steel	64 oz/in	7.00 N/cm	PSTC 101 Mod

Product Safety

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be noncytotoxic, non-sensitizing, and non-irritating.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60°F-80°F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY assumes any legal responsibility.



25 Forge Pkwy Franklin, MA 02038 1-800-343-7875 Fax 1-800-328-4822 www.berryglobal.com